


SCIENCE NEEDS DREAMERS

## CrestOptics X-Light Spinning Disk Scanner

CrestOptics X-Light marks an evolution in Nipkow spinning disk confocal scanner. Coupling with modern high intensity multi-wavelength solid state illumination sources, led sources and high resolution camera, X-Light enables high resolution confocal imaging in an easy to use and cost effective optical package that fits on your existing microscope.



High speed multi-point confocal scanning, combined with high speed sCMOS cameras, minimizes photobleaching and allows real-time imaging and recording at 500 fps or more. Moreover in combination with high sensitive EMCCD cameras allows low signal samples imaging. A long life light source coupled to the instrument allows for full visible spectrum (350nm - 750nm) confocal imaging (and NIR spectrum in some variants). Automation of internal multi-position excitation, dichroic and emission filter wheels allows fast multi-dimensional (X, Y, Z, t) imaging of up to five or more fluorescent probes in the same sample.




## Proprietary Spinning Disk Design

### Proprietary design high transmission low cross talk multi-spiral pinhole pattern


CrestOptics X-Light is a 15,000 RPM spinning disk with choice of single or double patterns spinning disk allowing users to select the appropriate pinhole size to match the objective lens numerical aperture or experiment protocol used.

### Long-term stability

X-Light systems show a long-term (years) stability with no need of further alignments after the first installation.



Convallaria stem, O. Garner Bergman-Labora Sweden.



## Microscope In-coupling

### Compatibility with a large variety of inverted and upright microscopes and full automation.

 It can be installed on all inverted and upright fluorescent microscopes model at 100% light port through C-mount adaptor (or F-mount adaptor in case of Large field of views variants). Thanks to the full motorized and software controllable X-Light filter wheels, even manual microscope are compatible for multi-color confocal imaging. Z-stack collection is done integrating manual microscope with a z-focus step motor.



• A motorized bypass mode between widefield and confocal mode, pinhole size selector, dichroic mirror wheel and emission filter wheel make X-Light easy to automate for a variety of experiments. Easy emission filter and dichroic filter exchange with the provided tools.

## Light Source In-coupling

Utilizing a solid state laser or LED based illumination source, the X-Light supports a wide array of excitation wavelengths with very long lifetimes. Illumination input to the X-Light is multimode fiber through an excitation gimbal mount for easy alignment on custom microscope setup and for the best S/N, producing flat, even illumination in the field-of-view. These light sources are fast switching for dynamic multichannel experiments, while still maintaining a very good investment to own and maintain.

## Camera In-coupling

## Optical relay lenses support up to 22mm FOV in standard configuration and 25mm FOV in large field of view configuration

 A wide selection of high-end cooled EM CCD cameras as well as large format sCMOS camera with a combination of 12-16 bit information, fast readout, high quantum efficiencies and small pixel size produces image at a high resolution and high signal-to-noise ratio.



### Optional image output splitter third party

- Optional TwinCam (splitter for 2 cameras dual imaging) or MultiCam (splitter for up to 4 cameras 4 channels imaging).
- Optional OPTOSPLIT II I or 2 images on a single camera: enabling a single camera to record images simultaneously at two different optical wavelengths or other differentiated state.



## Software Driver

### Full automation into a complete workstation

• Easy to operate: controls all motorized parts, automate image capture sequence.

### Choices of different package

 MetaMorph, NIS Elements or Micro-Manager plug-in for the X-Light makes it a very easy-to-learn and easy-to-use device for imaging. The X-Light is compatible with multidimensional acquisition, large image XY stitching applications, well plate scanning (HCA) and can be integrated on multimode systems to share a microscope stand with point scanners or other imaging systems.



First Generation Spinning Disk Confocal System LED and LASER Sources compatible.

- o Selected/high quality low autofluorescence optical components and anti-reflection coating treatments
- o Full compatibility with all inverted and upright microscopes
- o Full compatibility with third party components such as:
  - o External automated filter wheels
  - Optical Splitters





### Variants:

- X-Light VI "VIS": Standard system operating in the visible range 385nm-750nm
- X-Light VI "VIS-NIR": Standard system optimized in the range 385nm-1000nm
- X-Light VI "UV-VIS": System optimized in the UV-VIS range 350nm-750nm
- X-Light VI "NIR": System optimized for the near infrared fluorescent probes operating in the range 700nm-1600nm
- X-Light VI "20krpm": Standard system with increased disk speed up to 20000rpm



## CrestOptics X-Light V1 Specification

Light source —	Multichannel LED and Laser source compatibility Illumination source supports hardware triggering for fast multichannel experiment				
Supported input fiber	0.39NA multimode 1.5mnn fiber with SMA adapter; excitation gimbal mount for easy alignment on custom microscope setup and for best S/N				
Acquisition modes	- Widefield microscopy - Confocal microscopy - Bright-field microscopy - Phase-contrast microscopy				
	Double pattern disk 10mm x 10mm FOV each pattern for CCD	40 µm			
Pattern configuration		70 µm			
	Single pattern disk 22mm FOV for large format sCMOS	60 µm			
<u> </u>	Option: custom pattern available on request				
Disk speed	15,000 RPM standard (optional 20,000 RPM)				
Confocal resolution	With 60x NA 1.42 oil immersion Objective: <800nm				
Laser clean up filter	Option 3-position manual slider				
Dichroic wheel	Motorized 3-position (motorized 5-position on request)				
Dichroic size	Ultraflat 25.5 mm x 36 mm x 1 mm				
Emission filter wheel	Manual 4-position wheel standard (motorized 8-position on reques	t)			
Emission filter size	25mm diameter, up to 5mm thickness	7012,200000000 7012,20000000000000000000000000000000000			
Supported microscope	Upright and inverted microscope models from all brands with 100% c-mount	output port			
	High sensitivity CCD, sCMOS and EMCCD	72.5000000			
Camera	Easy camera focus adjusting internal optics without moving camera, no further disk and camera alignment needed				
Software control	Micro-Manager, MetaMorph, NIS Elements				



### Second Generation Spinning Disk Confocal System

- · Selected/high quality low autofluorescence optical components and anti-reflection coating treatments
- Full compatibility with all inverted and upright microscopes
- Full compatibility with third party components such as:
  - o External automated filter wheels
  - Optical Splitters

### Unique features compared to previous X-Light generation:

### • Plug-in Spinning Disk Box (possibility to change spinning disk in few seconds)

- Fast and easy opto-mechanical system for the excitation alignment
- New optical system design: diffraction limited over 22mm field of view
- New optical system design: diffraction limited over a large variety of objectives, from high magnification/high NA to low magnification/low NA objectives
- Full compatibility with Crest VCS (Video Confocal Super-resolution) module for 3D resolution enhancement: three microscopy systems at once





### Variants:

- X-Light V2 "VIS": Standard system operating in the visible range 385nm-750nm
- X-Light V2 LFOV 25mm FOV for Nikon Ti2 inverted microscope

## CrestOptics X-Light V2 Specification

Light source	Multichannel LED and Laser source compatibility			
	Illumination source supports hardware triggering for fast multichannel experiment			
upported input fiber	0.39NA multimode 1.5mm fiber with SMA adapter; excitation gimbal mount for easy alignment on custom microscope setup and for best S/N			
Acquisition modes	- Widefield microscopy - Confocal microscopy - Bright-field microscopy - Phase-contrast microscopy			
	Double pattern disk at 10mm x 10mm FOV each pattern for CCD			
Disk pinhole size vs camera size	Single pattern disk at 22mm FOV for large format sCMOS			
	Option: custom pattern available on request			
Disk speed	I 5,000 RPM standard			
Confocal resolution	60x NA 1.42 oil immersion Objective: <650nm			
Dichroic wheel	Motorized 5-position standard			
Dichroic size	Ultraflat 25.5 mm x 36 mm x 1 mm or 25.5mm x 36 mm x 2mm			
mission filter wheel	Motorized 8-position wheel standard			
Emission filter size	25mm diameter, up to 5mm thickness			
Excitation wheel	Manual 4 positions standard			
Excitation wheel	Motorized 8 positions (on request)			
ipported microscope	Upright and inverted microscope models from all brands with 100% c-mount output port			
	High sensitivity CCD, sCMOS and EMCCD			
Camera	Easy camera focus adjusting internal optics without moving camera, no further disk and camera alignment needed			
Software control	Micro-Manager, MetaMorph, NIS Elements			

## Optical Performances – V2 Diffraction Limited Performances

Reference wavelength: 587nm				
Magnification	NA	Theoretical Airy Disk Radius (um)	Actual Spot Size Radius (um)	Spot Size View*
100	1.45	24.8	4.1	
60	1.4	15.6	5.9	
20	0.75	9.6	7.6	


\* Black circle: theoretical Airy Disk




### Super-resolution imaging for any fluorescent stained sample

VCS (VideoConfocal super-resolution) module is a new system developed as an add-on for the X-light V2.

Using an approach based on super-resolution patents which exploit the 2-dimensional scan of the sample with a pinhole pattern and several algorithms application, CrestOptics-VCS enhance 3D resolution of the standard light microscope.





## Super-Resolution

## Three microscopy techniques in one solution

The full system X-Light V2 with VCS is thought to work in three configurations:

- Widefield mode
- Confocal mode
- Super-resolution mode

VCS technique includes the following phases:

- Specimen illumination with a multipoint beam obtained filtering the excitation light with a mask

- Multipoint beam scan along orthogonal directions (u,v) parallel to the illumination plane

- Acquisition of detected light coming from the specimen for each position of the illumination grid.

- Each acquired image is described by the distribution intensity for each pixel (x,y)  $I_{u,v}(x,y)$ 


- A final super-resolved image is calculated from the N raw images using super-resolution algorithms

## CrestOptics VCS Specifications

VCS resolution*	Lateral ~ 120 nm Axial ~ 300nm		
Axial Range	Up to 110 $\mu$ m		
Imaging	- 2D SIM - 3D SIM - Bypass mode for Confocal and Widefield acquisition		
Multi-Color Imaging	Motorized Focus Lens correction for multicolor imaging		
VCS Source	LED and Laser source compatible with multimode 1.5mm fiber, 0.39NA		
VCS Field of View	Up to 22mm FOV		
VCS objectives	Plan Apo corrected objectives, High Magnification, High NA compatibility: - 60x - 100x		
Operational conditions	No specific conditions required in temperature and sample preparation. Anti-vibration setup required		
VCS detector	sCMOS detectors high QE>80%, low readout noise required		
Software compatibility	Metamorph, Nis-Elements		
GPU video card/CUDA specification: - CUDA computational capability 5.0/5.2 - GPU Memory: 4 GB or more			

\* Depending on optics and super-resolution algorithms used

	A DESCRIPTION OF A	A STATE AND A STATE AND A STATE
Widefield	Confocal	vcs



Using three different microscopy acquisition modes offered by CrestOptics confocal and superresolution system, the resolution improvement of 3D reconstruction of  $1\mu$ m extended objects is shown.

# -FOV

Top Performance




Second Generation Spinning Disk Confocal System with enhanced field of view capability

- Unique confocal microscope to cover 25mm field of view
- Tailored to Nikon Eclipse Ti-2 inverted microscope
- · Selected high quality low autofluorescence optical components and anti-reflection coating treatments

Unique features with compared to previous X-Light generation:

- Plug-in Spinning Disk Box (possibility to change among multiple pinhole patterns in few seconds)
- · Fast and easy opto-mechanical system for the excitation alignment
- New optical system design: diffraction limited over 25mm field of view
- New optical system design: diffraction limited over a large variety of objectives, from high magnification/high NA to low magnification/low NA objectives
- Thanks to the extended field of view capability it acquires very large tissue areas with a strongly reduced number
  of acquisitions.



## CrestOptics X-Light V2 L-FOV Specification

Light source	Multichannel LED		
	Illumination source supports hardware triggering for fast multichannel experiment		
Supported input fiber	3mm LLG adapter; excitation gimbal mount for easy alignment on custom microscope setup and for best S/N		
Acquisition modes	- Widefield microscopy - Confocal microscopy - Bright-field microscopy - Phase-contrast microscopy		
Disk pinhole size	Single pattern disk at 25mm FOV		
vs camera size	Option: custom pattern available on request		
Disk speed	15,000 RPM standard		
Confocal resolution	60x NA 1.42 oil immersion Objective: <650nm		
Dichroic wheel	Motorized 5-position standard		
Dichroic size	Ultraflat 25.5 mm x 36 mm x 1 mm or 25.5mm x 36 mm x 2mm		
Emission filter wheel	Motorized 8-position wheel standard		
Emission filter size	25mm diameter, up to 5mm thickness		
	Manual 4 positions standard		
Excitation wheel	Motorized 8 positions (on request)		
Supported microscope	Nikon Eclipse Ti-2 inverted microscope		
	Large format cameras with F-mount		
Camera	Easy camera focus adjusting internal optics without moving camera, no further disk and camera alignment needed		
Software control	Micro-Manager, MetaMorph, NIS Elements		

## Confocal Imager

### **Combined simultaneous Photo-stimulation and Confocal Imaging**

X-Light Pro scanner system combines the disk scanner with a dual mirror galvo scanner for single point laser spot photo-stimulation. Control of the two independent scanning mirrors enables simultaneous stimulation and imaging to capture real-time reactions during stimulation at no delay time. Unsurpassed optical performance with selected high quality low autofluorescent optical components and anti-reflection coating treatments.

- FRAP Fluorescent Recovery after Photobleaching by photobleaching a small region of a sample in order to monitor the diffusion rate of fluorescent labeled molecules back into the photobleached region
- Uncaging
- Photo conversion protein
- · Optogenetic: receptor activation, protein recruitment, ion channel control
- Damage and repair






### **FRAP** features:

- Laser pulse time per scanning pixel > = 40us -Scanning bitmap ROI: up to 8 ROIs at the same time
- Scanning geomtrical ROI: up to 100 ROIs at the same time
- Autocalibration: 0 360 degree camera rotation -Software driver: MetaMorph
- Output signal power laser: 0 5 Vdc Angular Galvo mirror scaling: 10% 100%
- TTL laser pulse polarity inversion via serial command


## CrestOptics X-Light PRO for Confocal Imager



### Highly selective photobleaching scanning region

### Two types of ROI:

- Geometrical ROI: lines, squares, rectangles, points, multiline
- Bitmap ROI: Circles, ellipses, closed free-hand curves



## CrestOptics X-Light PRO for Confocal FRAP Specification

Software control	MetaMorph
	Dual fast galvo mirror scanning, 1.5ms sweep time, 256X256 or 512X512 points scanning resolution
	Single mode laser system (380 -780nm), laser pulse time per scan >=40us
FRAP features	Diffraction limit laser spot at sample <900nm diameter @470nm, scan up to 100 ROIs
	Software includes line, square, rectangle, point, multiline, circle, ellipse, closed free hand curves and areas
	Onboard firmware & fast RAM for automated operation

# **KEICHTNR**

•OP

First Generation Spinning Disk Confocal System working in the Near-Infrared spectrum

- Unique confocal microscope working in the range 700nm-1600nm
  - Excitation range: 700nm 1000nm
  - · Emission range: 800nm 1600nm

Unique features with compared to the standard X-Light generation:

• New optical system design: diffraction limited over 18mm field of view in the imaging range 900nm-1600nm

XILIGHTAN

• New optical system design: highly optimized AR coatings in the imaging range 800nm-1600nm

## CrestOptics X-Light NIR Specification


5


Light source	Laser			
Supported input fiber	0.39NA multimode 1.5mm fiber with SMA adapter; excitation gimbal mount for easy alignment			
Acquisition modes	- Widefield microscopy - Confocal microscopy - Bright-field microscopy - Phase-contrast microscopy			
	Double pattern disk 10mm x 10mm FOV each pattern 70µm			
Pattern configuration and Field of View	Single pattern disk 18mm FOV for large format      6	60µm		
	Option: custom pattern available on request			
Disk speed	15,000 RPM standard (optional 20000 RPM)			
Lateral resolution	25%-to-30% improvement with respect to widefield in the NIR region			
Laser clean-up filter	3-positions manual slider			
Dichroic wheel	Motorized 5-position standard			
Dichroic size	Ultraflat 25.5mm x 36mm x 1mm			
Emission filter wheel	Motorized 8-position			
Emission filter size	25mm diameter, up to 5mm thickness			
	Infrared camera in the 900nm-1600nm imaging range			
Camera	Easy camera focus adjusting internal optics without moving camera, no further disk and camera alignment needed			
Software control	Micro-Manager, MetaMorph, NIS Elements			
		State of the second		


## **Application Techinques**

### Routine and advance uses, ideal for living cell prolonged time lapse fluorescent imaging

• The significant advantages in the spinning disk approach are the ability to monitor rapidly occurring events within living cells without compromising resolution, as well as the high frequency low intensity signals because of substantially reduced photobleaching and phototoxicity.








E - Cell mitosis: Anti-tubulin antibody (green) and propidium iodide (red)


Investigation for Multicolor Colocalization, 3D optical sectioning and reconstruction and multi-dimensional imaging (X, Y, Z, t), especially ideal for live cell confocal imaging in much reduced photo-toxicity and photoleaching.

• Specimen fluorescent staining: fluorescent dye staining, immuno-fluorescence assay, fluorescent in-situ hybridization, fluorescent protein, microinjection of fluorescent labeled actin probes into living cells.



3D reconstruction skin sample, pan-neuronal marker, protein gene product 9.5 (Cy3), basement membrane marker, type IV collagen (CY2), endothelial cells are stained by Cy5-Ulex europeausagglutinin type I.


Dr. William R. Kennedy and Gwen Wendelschafer-Crabb, University of Minnesota, Minneapolis, MN



Time lapse recordings at 50 frames per second of calcium sparks in muscle cells loaded with calcium indicator dye Fluo-4 *Molecular Probes, Eugene, Oregon.* 


### Application areas

- Cell Biology & Plant Biology: apoptosis, autophagy, cell cycle, cell metabolism, cell tracking and tracing, cytotoxicity, oxidative stress detection, phagocytosis, endocytosis, receptor internalization, cell signaling, and communication, cell motility, cellular compartmentalization, protein synthesis and degradation, cellular and biophysical regulations.
- Cell and System Dynamics: structures and organs, e.g. blood vessels, neurons and processes such as angiogenesis and immune responses to vessel lesions
- Embryology & Developmental biology: C. elegans, Drosophila and Zebrafish growth and signal mechanisms
- Cancer research
- Clinical & translational medicine research
- Cardio and neuro sciences
- Calcium imaging, other ion measurements & membrane potential
- Yeast and bacteria studies
- Stem cell research and 3D cultures
- Synthetic biology: biofuels, vaccine & antibody production, plant sciences, industrial enzymes, biobased chemicals Pharma, biopharma & CRO




Muntjac cells with MouseAnti-OxPhos Complex V Inhibitor Protein, Alexa Fluor<sup>®</sup> 555 GoatAnti-Mouse IgG, Alexa Fluor<sup>®</sup> 488 Phalloidin and TO-PRO<sup>®</sup>-3

Widefield



Confocal





				Notes
X X Y X				
			20000000000000000000000000000000000000	
		00000000000000000000000000000000000000	10000000000000000000000000000000000000	
			10000000000000000000000000000000000000	20000000000000000000000000000000000000
		20000000000000000000000000000000000000	2000202000 20000000000	10000000000000000000000000000000000000
		7000000	50000000000000000000000000000000000000	00000000000000000000000000000000000000
				10000000000000000000000000000000000000
			70000000000000000000000000000000000000	00000000000000000000000000000000000000
				10000000000000000000000000000000000000
			,	700000000 70000000
				00000000000000000000000000000000000000
				0000000
			3000000	1000000
				0000000

## WorldWide Distribution

### EASTERN AMERICA, CENTRAL AMERICA, WESTERN AMERICA, SOUTH

Chris Baumann - Sales Manager

**89 North** Toll free 1.877.417.8313 Office + 1.802.881.0302 ext. 101 Mobile + 1.520.495.9661 Fax + 1.802.881.0308 E-mail: cbaumann@89north.com www.89north.com

Michael A. Davis - Bioscience System Product Manager Nikon Instruments Inc. Office: 631-547-8544 Fax: 631-547-4033 Mobile: 214-282-3837 E-mail: mdavis@nikon.net www.nikoninstruments.com

### **JAPAN** and **SINGAPORE**

Molecular Devices Tokyo head office Office: +81-3-6362-5260 (main), 5264(UIC team) Osaka Office Office: +81-6-7174-8831 E-mail: MetaMorph.Japan@moldev.com - Yasutomo.Kubota@moldev.com www.moleculardevices.com

#### **SOUTH KOREA**

Molecular Devices Tokyo head office Office: +81-3-6362-5260 (main), 5264(UIC team) Osaka Office Office: +81-6-7174-8831 E-mail: MetaMorph.Japan@moldev.com - Yasutomo.Kubota@moldev.com www.moleculardevices.com

YB.Yoon, Sales Manager VisionLink Corporation Office : +82-70-8752-1800 Fax : +82-31-349-3440 Email: info@visionlink.co.kr

### **AUSTRALIA and NEW ZEALAND**

Con Sapounas SciTech Pty Ltd - Imaging Spectroscopy Nanotechnology Melbourne office Office: +61 3 9480 4999 Mobile: +61 4 0044 0499 Sidney office Office: +61 2 9705 8059 (Sydney) Mobile: +61 3 9416 9959 E-mail: con@scitech.com.au www.scitech.com.au

### INDIA

Ramesh R - Chairman & Managing Director Medi Analytika Office: +91-44-24460988, 24463084, 24463085 Fax: +91-44-24463931 E-mail: info@medianalitika.com www.medianalitika.com

### Sanjeev Kumar

Towa Optics (I) Pvt. Ltd. Office : 011-26831936/4304 E-mail: sanjeev@towaoptics.com www.towaoptics.com

#### **UNITED KINGDOM**

Jeremy Graham - Managing Director Cairn Research (UK) Ltd. Toll Free: 0845 330 1267 (UK Only) Office: +44 (0) 1795 590140 E-mail: j.graham@cairn-research.co.uk www.cairn-research.co.uk

### ITALY

Alessandro Rossi - Biology Application Specialist Crisel Instruments s.r.l. Office: +39 06 35402933 E-mail: rossi@crisel-instruments.it www.crisel-instruments.it

### CHINA, TAIWAN, HONK KONG and MACAU

Harmony Group (H.K.) Limited Office: +852-96347908 Fax: +852-30139323 E-mail: info@harmonygrouphk.com www.harmonygrouphk.com

### EUROPE, ISRAEL, RUSSIA, AFRICA

Nikon Instruments Europe BV Laura Schneider - European Product Support Manager E-mail: Laura.Schneider@nikon.com Christelle Rosazza - European Product Specialist E-mail: Christelle.Rosazza@nikon.com Phone: +31207099500 www.nikoninstruments.eu

For further information please visit our website:

## www.crestopt.com

**CrestOptics s.p.a.** - Via Mattia Battistini, 184/D - 0016 - Roma - ITALY Ph. 0039.06.61660508 - Fax 0039.06.61662738 E-mail: info@crestopt.com