


# TRITOM

Multimodal whole body in-vivo Imaging Platform




Co-registered photoacoustic/ fluorescence 3D tomography

Molecular & functional imaging with high-resolution anatomical registration

Optimized for small murine research models




## **Co-registered Multimodality 3D Imaging**



True 3D anatomy with sub-millimeter spatial resolution:


- 160  $\mu$ m x 160  $\mu$ m in transverse planes
- 160 μm x 470 μm in sagittal and coronal planes




TriTom composite 890 nm + 532 nm imaging of a life mouse vasculature


#### **Enhance imaging benefits, eliminate shortcomings**

 TriTom<sup>TM</sup> utilizes simultaneous spectrallyselective optical and photoacoustic excitation and detection to create unparallel volumetric assessment of live organisms, organs, and tissues.






Molecular sensitivity of photoacoustics is enhanced by fluorescence




Fluorescence of a contrast agent (780/800 nm excitation/emission) + photoacoustic skin vasculature

# **TriTom Applications**

#### Triple Analysis – anatomical, functional, and molecular




#### **Cancer Research**




Composite skin & deep tissue 3D image. 532 nm skin excitation. 890 nm deep tissue excitation. Tumor size 10.6 x 4.7 x 11.6 mm3.

#### **Neuroscience**




7 mm MIP axial slab of mouse's brain near the cerebellum/medulla. (1) confluence of sinus and (2) cerebral artery are marked.

#### **Functional Imaging**




Lower abdomen of a healthy mouse. Red – deep blood vessels, 890 nm excitation scan. Yellow – superficial vasculature (skin), 532 nm excitation scan. 1 – right subiliac lymph node. 2 – injection site (right mammary fat pad)

#### **Stem Cell Research**




3D photoacoustic image of the spinal cord with injected cells denoted by arrows.

#### **Developmental Biology**



Pregnant mouse, gestational day 12 shows a composite skin / deep tissue photoacoustic image.

#### **Anatomical Details of Internal Organs**



PAT maximum intensity projection of mouse liver. (1) left lobe of liver; (2) right lobe of liver; (3) hepatic vein branches.


#### **Development of Contrast Agents**



Crossectional (axial) views of 0.8 mm inner diameter of the microcuvettes with IRDye800CW samples.


#### **Accessories**

# Cylindrical cuvette holder (development of contrast agents)



- High throughput:
- ⇒ Interrogate up to 10 samples per scan
- ⇒ Preparation time < 5 min
- Save your valuable samples!
- ⇒ Tiny 50 µL sample volumes
- Convenient & repeatable procedure:
- ⇒ Radially oriented slots consistently hold samples in place, allow quick setup and removal of cuvettes
- (1) Four  $\emptyset$  0.8 mm cylindrical cuvettes with samples
- (2) Plastic support rods
- (3) Radial slots for quick setup and removal of cuvettes
- (4) Silicon sealant applied to the ends of the cuvettes
- (5) Magnetic connects
- (6) Port for administering liquid scattering background

# Mouse holder (in vivo imaging)




- · Reliable gas anesthesia
- ⇒ Delivered through the hollow shaft
- Convenient & repeatable procedure:
- ⇒ Mouse's front and hind legs and the head are fixed in a consistent position with minimal stress
- ⇒ Preparation time < 5 min
- (1) Hollow shaft to fix the mouse holder inside TriTom while reliably delivering anesthesia gas
- (2) Cushioned attachment slots for mouse's front legs
- (3) Adjustable support block for mouse's hind legs with cushioned attachment slots
- (4) Bite bar
- (5) Plastic support rods

#### Additional options include:

Phantom kits

- Animal heating pad
- Animal anesthesia system
- Animal holder

### **Patent Protected Technology**



#### Perfect co-registration of imaging modalities

Both PA & FL are initiated simultaneously using pulsed high-efficiency laser light

Both PA & FL are enabled for molecular imaging by the tuning excitation wavelength during the scan. Tuning is necessary for multi-spectral imaging.

#### **Excitation Source**

Excitation source: Tunable OPO laser

Wavelength tuning range covering <u>all popular</u> <u>visible, NIR I and NIR II fluorophores and</u> <u>nanoparticles:</u>


- ⇒ 660 1064 nm (standard)
- ⇒ 460 659 nm (extended visible excitation)
- $\Rightarrow$  1065 1300 nm (extended NIR II excitation)

Up to 250 mJ peak energy @ 700 nm

10 Hz or 20 Hz pulse repetition frequency

Excitation linewidth < 0.5 nm (equivalent to 1,280 excitation filters!)

Integrated energy meter (quantitative imaging)





Photoacoustic and fluorescence excitation through the same ports (4)


### **Software**

The **TriTom**<sup>TM</sup> imaging software suite provides a streamlined experience from simultaneous acquisition of optical and photoacoustic data to final reconstruction and co-registration of volumes.




- Iconographic design of the main control screen
- Modular user-centric graphic Interface
- Optical frames
- Scan metadata
- Open data formats





TriTom image & signal monitoring interface

3D slicer — visualization and image processing software



# **Technical Specifications**

Navelength Range [nm]	460-1300	Single output, continuous tuning
Pulse Width [ns]	3 - 5	FWHM
Linewidth [cm-1]	< 10	
Repetition Rate [Hz]	10/20	
Pulse Energy [mJ]	> 120 mJ @ 700 nm	
before fiber bundle transmission	> 10 mJ @ 500 nm	
Pulse Energy Fluctuations [%]	< 2	StDev
Energy meter	Real-time pulse energy measurements stored in data header	
Fast wavelength switching	[470, 659] nm, [660, 1300] nm	
Dimensions [in (cm)]	26.5 (67.2) x 17.1 (43.4) x 34.9 (88.7)	
Weight [lbs (kg)]	150 (68)	
Input Power	208 or 240 V 10 A 50 / 60 Hz	

Fiber Bundle		
Input / Output Configuration	1:4 Circular/Linear	Over 2500 individually randomized fibers
Axial Excitation Spot Size [mm]	30 (50, 80)	standard (optional)
Length [m]	2	

Photoacoustic Transducer Array		
Configuration	Curved	Cylindrical focusing
Number of Elements	96	No dead elements
Element Size [mm^2]	1.3 x 1.3	Measured along centerline
Element Pitch [mm]	1.4	
Active angular aperture [deg]	0.1	
Radius of curvature [mm]	118	
Central Frequency [MHz]	6 ± 10%	T/R measurements, Optimized sensitivity in receive mode
Bandwidth @ -6 dB [%]	≥ 55	T/R measurements
Acoustical Matching	water	1.5 MRayl
Transducer Material	PEEK	Housing is connected to ground to prevent RF noise pickup
Utilization	Continuous immersion under 0.5 m of water between 10 to 40°C	
Shielding	Metalized inside, partial metalized outside (e.g. flash of gold)	

Fluorescence Emission Filters		
Optical Filter Wheel	Motorized	USB 2.0 PC connection
Clear Aperture [mm]	25	
Standard Optical Filter Emis-	11 filters covering emission range between 483 nm and 850 nm, 1 blocked, 1 open, and 1	
sion Ranges [nm]	custom slots	

Fluorescence Detector		
Detector Type	sCMOS	High sensitivity cooled scientific camera
Bit Depth	16-bit	
Sensor Size [in]	1.2	
Number of Pixels	2048 x 2048	
Pixel Size [μm]	6.5 x 6.5	
Quantum Efficiency [%]	20 - 95	200 - 1000 nm
Readout Noise [e-]	1.2	Low readout noise for high frame rate applications
Dark Current [e-]	0.03	For 100 ms or shorter exposures
Dark Signal Nonuniformity [e-]	0.2	
Max Frame Rate [fps]	35	Full resolution
Cooling [°C]	-15	Peltier cooling
PC Connection	USB 3.0	

GION™ ADC Data Acquisition	Unit	
Channels	256	
Programmable Gain [dB]	46 - 91	Measured using oscilloscope with 50 $\Omega$ input
Analog Bandwidth	40 kHz - 25 MHz	
Resolution	12-bit	
Sampling Rate [MSPS]	40	
Min Input Impedance $[k\Omega]$	1	
Max Frame Rate [Hz]	200	Up to 400 Hz using 128-channels only
Points / Frame / Channel	4096	
Trigger Connections	2x SMA Electrical / 2x Optical	

Form Factor	Desktop	Mini ITX case with handle	
Processor	Core i5-10600	6 Cores / 12 Threads 3.3 - 4.8 GHz 12 MB 65 W	
Graphics	Geforce RTX 2070 Super	2560 CUDA cores 1605 - 1770 MHz 8 GB GDDR6 215 W	
Memory	16 GB	DDR4 3200 MHz	
Storage	2 TB	M.2 NVMe PCIe 3.0 SSD	
Imaging Monitors	Two 27" Monitor, resolution	Two 27" Monitor, resolution 2048x1080	
Operating System	Microsoft Windows 10 Pro	Microsoft Windows 10 Pro 64-bit	
Dimensions [in (cm)]	14.8 (37.5) x 7.83 (19.9) x 1	14.8 (37.5) x 7.83 (19.9) x 12.3 (31.2)	
Weight [lbs (kg)]	20 (9.1)		
Imaging Software	Advanced TriTom imaging suite (multimodality 3D imaging, molecular spectral imaging, molecular unmixing); 3D slicer visualization & image analysis		



# **About PhotoSound ®**

PhotoSound Technologies, Inc. (Houston, Texas USA) develops new imaging products and technologies for life sciences. A 3D imaging platform for in vivo preclinical research and drug discovery (TriTom™) is implemented on patented PhotoAcoustic Fluorescent Tomography (PAFT) technology, which utilizes simultaneous spectrally-selective optical and photoacoustic excitation and detection to create unparallel volumetric assessment of live organisms, organs, and tissues. A MoleculUS™ system is developed for clinical research that can benefit from co-registered ultrasound and molecular photoacoustic imaging. We also offer a variety of OEM electronic components for multi-channel parallel data acquisition.

All PhotoSound technology solutions are designed and built by experts in biomedical imaging systems, photoacoustics, ultrasound, optics, electronics and tunable lasers. Our employees are committed to provide every customer with the highest quality products and services, short delivery times and competitive pricing. Visit us at www.photosound.com to learn more about our products and proprietary technologies.

Version TRITOM.003.0620 © 2021

All specifications are subject to change without notice.

**TriTom**<sup>TM</sup> is classified EAR99 and does not require an export license.

PhotoSound® is a registered trademark of PhotoSound Technologies, Inc.

**TriTom**<sup>TM</sup> is designated for pre-clinical research only and is not suitable for human clinical trials.

PhotoSound \*

PhotoSound Technologies, Inc. | Imaging and Data Acquisition Solutions
9511 Town Park Drive | Houston, TX 77036 USA