

Optolong Fluorescence Filters



Item No.	Product Name	Specs.	Size
30024	DAPI	EX BP365/50	φ 25 x 5mm
		EM BP445/50	φ 25 x 3.5mm
		DM 405	25.7 x 36 x 1.05mm
30035	DAPI-2	EX BP375/28	φ 25 x 5mm
		EM BP460/50	φ 25 x 3.5mm
		DM 415	25.7 x 36 x 1.05mm
30025	FITC	EX BP475/40	φ 25 x 5mm
		EM BP530/50	φ 25 x 3.5mm
		DM 500	25.7 x 36 x 1.05mm
30016	FITC-2	EX BP480/30	φ 25 x 5mm
		EM BP535/40	φ 25 x 3.5mm
		DM 505	25.7 x 36 x 1.05mm
30017	TRITC	EX BP545/20	φ 25 x 5mm
		EM BP605/60	φ 25 x 3.5mm
		DM 560	25.7 x 36 x 1.05mm
30027	TRITC-2	EX BP545/30	φ 25 x 5mm
		EM BP605/40	φ 25 x 3.5mm
		DM 575	25.7 x 36 x 1.05mm
30026	TexasRed	EX BP560/30	φ 25 x 5mm
		EM BP630/70	φ 25 x 3.5mm
		DM 585	25.7 x 36 x 1.05mm
30021	Cy5	EX BP630/40	φ 25 x 5mm
		EM BP695/40	φ 25 x 3.5mm
		DM 660	25.7 x 36 x 1.05mm
30031	FISHAqua	EX BP440/20	φ 25 x 5mm
		EM BP480/30	φ 25 x 3.5mm
		DM 460	25.7 x 36 x 1.05mm
30032	FISHGreen	EX BP490/20	φ 25 x 5mm
		EM BP530/25	φ 25 x 3.5mm
		DM 505	25.7 x 36 x 1.05mm
30033	FISHOrange	EX BP540/20	φ 25 x 5mm
		EM BP590/25	φ 25 x 3.5mm
		DM 560	25.7 x 36 x 1.05mm
30034	FISHRed	EX BP580/20	φ 25 x 5mm
		EM BP630/25	φ 25 x 3.5mm
		DM 605	25.7 x 36 x 1.05mm

Item No.	Product Name	Specs.	Size
31001	LF-DAPI(DAPI longpass)	EX BP365/50	ϕ 25 x 5mm
		EM LP415	ϕ 25 x 3.5mm
		DM 405	25.7 x 36 x 1.05mm
31003	LF-FITC(FITC longpass)	EX BP475/40	ϕ 25 x 5mm
		EM LP510	ϕ 25 x 3.5mm
		DM 500	25.7 x 36 x 1.05mm
31004	LF-TRITC(TRITC longpass)	EX BP540/20	ϕ 25 x 5mm
		EM LP590	ϕ 25 x 3.5mm
		DM 565	25.7 x 36 x 1.05mm
31005	LF-TexasRed(TexasRed longpass)	EX BP530/50	ϕ 25 x 5mm
		EM LP590	ϕ 25 x 3.5mm
		DM 575	25.7 x 36 x 1.05mm
32012	Dual-bandGreen/Orange	EX BP480-560	ϕ 25 x 5mm
		EM BP525-605	ϕ 25 x 3.5mm
		DM 500-580	25.7 x 36 x 1.05mm
32013	Dual-bandFITC/TexasRed	EX BP480-575	ϕ 25 x 5mm
		EM BP525-635	ϕ 25 x 3.5mm
		DM 500-600	25.7 x 36 x 1.05mm
33011	Triple-bandDAPI/Green/Orange	EX BP395-480-560	ϕ 25 x 5mm
		EM BP440-521-606	ϕ 25 x 3.5mm
		DM 410-500-580	25.7 x 36 x 1.05mm
33012	Triple-bandDAPI/FITC/TexasRed	EX BP405-495-573	ϕ 25 x 5mm
		EM BP460-533-630	ϕ 25 x 3.5mm
		DM 410-510-590	25.7 x 36 x 1.05mm

Item No.	Product Name	Specs.	Size
Olympus U-MF2		<p>Compatible with Olympus :</p> <p>AX Series - AX70</p> <p>BX2 Series - BX41, BX50, BX51, BX60, BX61, BX50/51WI, BX60/61WI</p> <p>IX2 Series - IX50, IX51, IX70, IX71, IX81</p>	
Olympus U-FF		<p>Compatible with Olympus :</p> <p>BX3 Series - BX53, BX63</p> <p>IX3 Series - IX53, IX73, IX83</p>	
Zeiss		<p>Compatible with ZEISS :</p> <p>Axio Imager, Axiostar plus, Axioskop40, Axioskop2(post-2001), Axioplan2i, Axio2ie, Axiovert200, Axiovert40, Axio Observer, Axio Examiner, AxioScope A1</p>	
Nikon		<p>Compatible with Nikon :</p> <p>TE2000, 50i, 55i, 80i, 90i, Eclipse Ti, Ni, Ci series</p>	
Leica		<p>Compatibe with Leica :</p> <p>DM2000, DM2500, DM3000, DI3000B, DM4000, DMI4000B, DM5000, DM5500, DM6000, DMI6000B, DM4, DM6</p>	