



μ Manager

For Compatibility: See Software Provider							
Illumination System	Recognised in µManager as:	Required µManager version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	1.4 or later	Individual Wavelength control	Individual Wavelength control	N/A	<u>Click for link to</u> <u>Micromanager</u> <u>website</u>	Windows XP onwards. 32-bit & 64-bit <u>Click for link to</u> <u>CoolLED driver</u>
pE-300 ^{white}	pE-300	1.4 or later	Individual Wavelength control	Global control	N/A	pE-300 & pE-4000	Windows 64-bit only
pE-300 ^{ultra}	pE-300	1.4 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A	device adapters are present in μManager <u>1.4.23</u> nightly builds as standard	<u>Click for link to pE-</u> <u>300.dll</u>
pE-4000	pE-4000	1.4 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control	The device adapters can also be downloaded and	Windows 64-bit only <u>Click for link to pE-</u> <u>4000.dll</u>
pE-800	pE-800	1.4	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control	added to the official <u>1.4.22</u> release only.	Windows 64-bit only <u>Click for link to pE-</u> <u>800.dll</u>
pE-340 ^{fura}	Not recognised	N/A	N/A	N/A	N/A	N/A	N/A
pE-6501/ pE-6501-8	NI100X	1.4 or later	N/A	N/A	N/A	<u>pE-6501 Guide</u> <u>pE-6501-8 Guide</u>	N/A